



Viruclear VF Series Virus Removal Filter





◀ Viruclear VF Series plastic housing cassette



Viruclear VF

Reliable & Effective Virus Removal Filter

Virus filtration is a robust virus clearance method in the production process of common biopharmaceutical products. Due to its simple and gentle operation, clear mechanism, and easy verification, it has been widely used in the downstream purification process of blood products such as coagulation factors, immunoglobulins, and recombinant proteins and antibody biological products expressed in animal cells to ensure virus safety.

The most widely used virus filtration product is the 20nm pore size parvovirus removal filter in the current market. Size exclusion is the main mechanism of virus removal. The target products with molecular sizes smaller than the pore size of the membrane can pass through the membrane, and the virus with sizes larger than the membrane pore size can be intercepted.

Based on its physical interception characteristics, whether the diameter of 80-120 nm mouse leukaemia virus is a retrovirus or the diameter of about 20 nm of parvovirus, can be reliable and efficient interception. Virus clearance meet regulatory requirements that log reduction values (LRV) \geq 4 log.

Cobetter Viruclear VF series virus removal filters features an asymmetric polyethersulfone (PES) membrane, greatly improve the fouling space, and provides strong interception and removal fouling for protein aggregates and potential viruses. The part near the water surface consists of a precision retention layer with about 20 nm and the pore size distribution is uniform, which can achieve robust retention of parvovirus and robust retention maintained during process interruption/depressurization. Cobetter Viruclear VF series virus removal filters also provide good scalability, while integrity testing using an easy-to-operate diffusion flow method.

Features

- Unique surface hydrophilic improvement of PES membrane , which gives high mass capacity and high yield of virus retentive membrane
- Asymmetric membrane structure, high porosity, and high flux
- Robust removal of parvovirus
- \geq 4.0 log removal of model virus
- Robust retention maintained during process interruption
- Chemically stable
- Easy to install, use and test

Typical Applications

- Monoclonal Antibody
- Bispecific Antibody (BsAb)
- Antibody-Drug Conjugate(ADC)
- Fusion Proteins
- Nanobody
- Recombinant Proteins


Viruclear™ VF

PES Virus Removal Filters for Antibody Processing

Cobetter Viruclear™ VF virus removal membrane utilizes a highly hydrophilic formulation, effectively enhancing the hydrophilicity of the PES membrane. The Viruclear™ VF filter enables robust linear process scalability and supports inline integrity testing. It is designed primarily for purification applications of biologics, including mammalian cell-expressed recombinant proteins, mAbs, ADCs, and high-purity/stability bispecific antibodies.

Excellent Product Adaptability

The virus removal performance of protein solutions is strongly correlated with the intrinsic properties of the solution. Product with high protein concentration, low purity, high hydrophobicity, or containing stabilizing agents typically pose greater challenges to virus filtration. The figure presents representative examples of the filtration performance of Viruclear™ VF when processing different customer protein solutions.



Viruclear™ VF achieves a mass loading of 5–16 kg/m² across various customer mAb products during a 3-hour filtration at a constant pressure of 2.0 bar.

This demonstrates that Viruclear™ VF offers high mass throughput, meeting the filtration loading requirements for most mAbs.

Robust Removal of Virus


Test Condition

Viruclear VF syringe filters performed MVM challenge experiments on multiple customer mAb batches at room temperature and 30 psi.

Filter	Viruclear VF syringe filter
Model virus	MVM
Feed solution	mAb protein solutions from different customers
Constant pressure	30 psi
Flux decay	Decay to 25% of initial flux or target loading achieved



Test results



Viruclear™ VF product possesses robust virus clearance capability, effectively ensuring stable virus retention exceeding regulations. Meanwhile, addressing worst-case conditions in production—pressure interruption, Viruclear™ VF still ensures robust virus retention

Viruclear™ VF Plus

Strength-Optimized PES Virus Removal Filter



Through optimized membrane fabrication, Cobetter introduces Viruclear™ VF Plus, featuring enhanced structural strength and improved manufacturing yields. For certain products, VF Plus delivers superior filtration performance.


Both Viruclear™ VF Plus and Viruclear™ VF are designed for virus filtration in biologics, including recombinant proteins, mAbs, ADCs, and high-purity bispecific antibodies. Users should select the most suitable model based on process development test results.

Excellent Product Adaptability

The virus removal performance of protein solutions is strongly correlated with the intrinsic properties of the solution.

Product with high protein concentration, low purity, high hydrophobicity, or containing stabilizing agents typically pose greater challenges to virus filtration.

The figure presents representative examples of the filtration performance of Viruclear™ VF Plus when processing different customer feed materials.



Viruclear™ VF Plus achieves a mass loading of 6–15 kg/m² across various customer mAb products during a 3-hour filtration at a constant pressure of 2.0 bar.

This demonstrates that Viruclear™ VF Plus offers high mass throughput, meeting the filtration loading requirements for most mAbs.

Robust Removal of Virus


Test Condition

Viruclear VF Plus syringe filters performed MVM challenge studies on multiple customer mAb batches at room temperature and 30 psi.

Filter	Viruclear VF Plus syringe filter
Model virus	MVM
Feed solution	mAb protein solutions from different customers
Constant pressure	30 psi
Flux decay	Decay to 25% of initial flux or target loading achieved



Test results



Viruclear™ VF Plus product possesses robust virus clearance capability, effectively ensuring stable virus retention exceeding regulations. Meanwhile, addressing worst-case conditions in production—pressure interruption, Viruclear™ VF still ensures robust virus retention


Viruclear VF Plus X

Enhanced Capacity PES Virus Removal Filter

Viruclear™ VF Plus X further enhances membrane dirt holding capacity. It has better adaptability to complex products and higher mass or volume capacity at the expense of a certain water flux, especially suitable for biological products with high protein aggregates content and poor stability, such as bevacizumab, unstable bispecific antibodies and fusion proteins. Viruclear VF Plus X is a supplement and alternative to Viruclear VF Plus in the fields of recombinant protein and antibody applications, and can take care of both filtration applications for some biological enzymes and low concentration blood products.

Enhanced Fouling Resistance

The comparative experiments on the filtration performance of Viruclear™ VF Plus and Viruclear™ VF Plus X were conducted under room temperature and constant-pressure filtration conditions of 30 psi, with parallel experiments filtering the same 0.4 g/L intravenous immunoglobulin (IVIG) solution used as the loading evaluation solution. The comparison of filtration flux decay trends can reflect the differences in filtration performance among the filters. The comparison of flux decay trends demonstrates the enhanced fouling resistance of Viruclear™ VF Plus X.



Robust Removal of Virus


Test Condition

Viruclear VF Plus X syringe filters performed MVM challenge studies on multiple customer mAb batches at room temperature and 30 psi.

Filter	Viruclear VF Plus X syringe filter
Model virus	MVM
Feed solution	mAb protein feedstreams from different customers
Constant pressure	30 psi
Flux decay	Decay to 25% of initial flux or target loading achieved



测试结果



Viruclear™ VF Plus X product possesses robust virus clearance capability, effectively ensuring stable virus retention exceeding regulations. Meanwhile, addressing worst-case conditions in production—pressure interruption, Viruclear™ VF Plus X still ensures robust virus retention


Virus Removal Filtration Process Development Service

Our technical engineers will work with customers to optimize the process parameters during the process development of virus removal filtration to obtain a robust, efficient and economical filtration operation.

Virus Removal Filterability Testing Device



Virus Removal Assembly Sketch Map



Specification

	Specification	Effective Filtration Area	Application
Viruclear™ VF Virus Removal Filter	DC series	2.8cm ²	For process development and virus clearance validation research
	Pilot scale	0.017/0.07/0.22m ²	≤200L
	Production scale	0.50/1.50m ²	>200L

Material of Products

	Specification	Membrane	Format	Shell	Effective Filtration Area	Assessries material
Viruclear™ VF Series	DC series	PES	Syringe filter	PVDF	2.8cm ²	—
	Pilot scale	PES	Single-use Plastic-housing Cassette	PVDF	0.017/ 0.07/ 0.22m ²	Silicone, PP
	Production scale	PES	Single-use Plastic-housing Cassette	PVDF	0.50 / 1.50m ²	Silicone, PP
Corevital™ SMDA Series	DS series	PVDF	Syringe filter	PP	0.72/3.4/4.5cm ²	—
	Pilot scale	PVDF	Single-use Plastic-housing Capsule	PP	180/420/660cm ² 0.12/0.16/0.28m ²	Silicone rubber
	Production scale	PVDF	Single-use Plastic-housing Capsule	PP	0.55/1.1/1.65m ²	Silicone rubber
Viruclear™ PNY Series	DS series	Nylon	Syringe filter	PP	3.4cm ²	—
	Pilot scale	Nylon	Single-use Plastic-housing Capsule	PP	0.025/0.12/0.3m ²	Silicone
	Production scale	Nylon	Single-use Plastic-housing Capsule	PP	0.6/1.2/1.8m ²	Silicone
Viruclear™ PCE Series	DS series	PES	Syringe filter	PVDF	2.8cm ²	PP
	Pilot scale	PES	Single-use Plastic-housing Cassette	PVDF	0.017/0.07/0.22m ²	PP
	Production scale	PES	Single-use Plastic-housing Cassette	PVDF	0.5/1.5m ²	PP
Viruclear™ PDT Series	DS series	Cellulose, diatomite, Nylon	Syringe filter	PP	4.5cm ²	—
	Pilot scale	Cellulose, diatomite, Nylon	Single-use Plastic-housing Cassette Single-use Plastic-housing Cassette(M02)	PP PP/ GF	0.027/0.15/0.4m ² 0.027/0.054m ²	Silicone, PP
	Production scale	Cellulose, diatomite, Nylon	Single-use Plastic-housing Cassette Single-use Plastic-housing Cassette(M02)	PP PP/ GF	0.92/1.1m ² 0.11/0.55/1.1m ²	Silicone, PP
Viruclear™ PDS Series	DS series	Silicon dioxide, Composite fiber, Nylon	Syringe filter	Glass fiber reinforced polypropylene	4.5cm ²	Silicone rubber
	Pilot scale	Silicon dioxide, Composite fiber, Nylon	Single-use Plastic-housing Cassette	Glass fiber reinforced polypropylene	0.027/0.054m ²	Silicone rubber
	Production scale	Silicon dioxide, Composite fiber, Nylon	Single-use Plastic-housing Cassette(M02)	Glass fiber reinforced polypropylene	0.11/0.55/1.1m ²	Silicone rubber

Ordering Information

Corevital™ SMDA Virus Removal Prefilter (PVDF Prefilter)

Bricap™ SFU Syringe Filters



Syringe filter

U13 Filtration area 0.72cm²

U25 Filtration area 3.4cm²

U33 Filtration area 4.5cm²



Connection Type

CP Luer lock female inlet and

Luer lock male outlet

SMDA

Filter Membrane Configuration

SMDA SMDA 0.2/0.1µm



Sterilization

G Gamma Compatible



Quantity/Package

1 1pcs / pk



Market

P Biopharmaceutical



Bricap™ C Series Capsule Filter



Housing type

Blank Normal

K Transparent



Syringe filter

C01 Filtration area 180cm²

C02 Filtration area 420cm²

C03 Filtration area 660cm²



Inlet/Outlet type

T 19mm(3/4")Sanitary flange

K 14mm(9/16")Sanitary flange

F 31mm(1 1/4")Sanitary flange

SMDA

Filter Membrane Configuration

SMDA SMDA 0.2/0.1µm



Sterilization

A Autoclavable only

G Gamma Compatible

S Sterile

C Sterile packaging with sterilizing breather bag



Quantity/Package

1 1pcs / pk



Product category

Blank Standard

Q

P

Market

P Biopharmaceutical



Bricap™ L Series Capsule Filter



Housing type

Blank Normal

K Transparent



Syringe filter

L02 Filtration area 0.12m²

L03 Filtration area 0.16m²

L05 Filtration area 0.28m²

L10 Filtration area 0.55m²

L20 Filtration area 1.1m²

L30 Filtration area 1.65m²



Inlet/Outlet type

S 38mm(1 1/2")Sanitary flange

T 19mm(3/4")Sanitary flange

K 14mm(9/16")Sanitary flange

SMDA

Filter Membrane Configuration

SMDA SMDA 0.2/0.1µm



Sterilization

A Autoclavable only

G Gamma Compatible

S Sterile

C Sterile packaging with sterilizing breather bag



Quantity/Package

1 1pcs / pk



Product category

Blank Standard

Q

P

Market

P Biopharmaceutical



Ordering Information

Viruclear™ PNY Virus Removal Prefilter (Nylon Filter)

V P

Application

VP Viruclear PNY

N Y

Material of Membrane

NY Nylon

Type

DS Syringe filter
(3.4cm²)

Quantity/Package

N9 Only for syringe filter,
9pcs / pk
N1 1pcs / pk

P

Market

P Biopharmaceutical



Viruclear™ PNY Virus Removal Prefilter (Nylon Filter)

V P

Application

VP Viruclear PNY

N Y

Material of Membrane

NY Nylon

Type

C02 250cm²
L02 0.12m²
L05 0.30m²

Type

TT 3/4" TC

Quantity/Package

N1 1pcs / pk

P

Market

P Biopharmaceutical



Viruclear™ PNY Virus Removal Prefilter (Nylon Filter)

V P

Application

VP Viruclear PNY

N Y

Material of Membrane

NY Nylon

Type

L10 0.60m²
L20 1.20m²
L30 1.80m²

Type

SS 1.5" TC

Quantity/Package

N1 1pcs / pk

P

Market

P Biopharmaceutical



Ordering Information

Viruclear™ PCE Virus Removal Prefilter with Plastic Shell

V P

Application

VP Viruclear PDF

C E

Material of Membrane

CE Cation Exchange
PES Membrane

□ □

Type

DC Syringe filter
(2.8cm²)

□ □


Quantity/Package

N9 9pcs / pk
N1 1pcs / pk

P

Market

P Biopharmaceutical



Viruclear™ PCE Virus Removal Prefilter with Plastic Shell

V P

Application

VP Viruclear PDF

C E

Material of Membrane

CE Cation Exchange
PES Membrane

□ □

Type

LA Lab(0.017m² & 0.007m²)
FL Flow(0.22m² & 0.5m² & 1.5m²)

□ □ □

Type

002 0.017m²
008 0.007m²
026 0.22m²
050 0.5m²
150 1.5m²

□ □


Quantity/Package

N1 1pcs / pk

P

Market

P Biopharmaceutical



Viruclear™ PDS Virus Removal Prefilter

V P

Application

VP Viruclear PDF

D S

Material of Membrane

DS Depth filter

□ □

Type

DS Syringe filter (4.5cm²)
27 M02(0.027m²)
54 M02(0.054m²)
M1 M02(0.11m²)
M5 M02(0.55m²)
MX M02(1.1m²)

□ □

Quantity/Package

N9 9pcs / pk
N1 1pcs / pk

P

Market

P Biopharmaceutical



Ordering Information

Viruclear™ PDT Virus Removal Prefilter (Depth Filter)

V P

Application

VP Viruclear PDT



D T

Material of Membrane

DT Depth filter

Type

DS Syringe filter
(4.5cm²)

Quantity/Package

N9 Only for syringe filter,
9pcs / pk
N1 1pcs / pk

P

Market

P Biopharmaceutical

Viruclear™ PDT Virus Removal Prefilter (Depth Filter)

V P

Application

VP Viruclear PDT



D T

Material of Membrane

DT Depth filter

Type

L4 L08TT(0.027m²)
SA CSCD(0.15m²)
SC CSCE(0.4m²)
27 M02(0.027m²)
54 M02(0.054m²)

Quantity/Package

N9 9pcs / pk
N1 1pcs / pk

P

Market

P Biopharmaceutical

Viruclear™ PDT Virus Removal Prefilter (Depth Filter)

V P

Application

VP Viruclear PDT



D T

Material of Membrane

DT Depth filter

Type

SB CSCB(0.92m²)
M1 M02 Cassette(0.11m²)
M5 M02 Cassette(0.55m²)
MX M02 Cassette(1.1m²)

Quantity/Package

N1 1pcs / pk

P

Market

P Biopharmaceutical

Ordering Information

Viruclear™ VF Virus Removal Syringe Filters with PVDF Shell

V F

Application

VF Virus filtration

E

Material of Membrane

E PES

S

Configuration Code

S Standard

□ □

Type

DC Syringe filter
(2.8cm²)

□ □

Quantity/Package

N9 Only for syringe filter,
9pcs / pk
(3 batches, 3 pcs per batch)
N1 1pcs / pk

P

Market

P Biopharmaceutical



Viruclear™ VF Removal Plastic-housing Cassette

V F

Application

VF Virus filtration

E

Material of Membrane

E PES

S

Configuration Code

S Standard

□ □

Type

LA Lab
(0.017/0.07m²)
FL Flow
(0.22/0.50/1.50m²)

□ □ □

Type

002 0.017m²
008 0.07m²
026 0.22m²
050 0.50m²
150 1.50m²

□ □

Quantity/
Package

N1 1pcs / pk

P

Market

P Biopharmaceutical



Ordering Information

Viruclear™ VF Plus Virus Removal Syringe Filters with PVDF Shell

V F

Application

VF Virus filtration

P

P Plus

E S

Material of Membrane

ES mPES,
hydrophilic modified
PES membrane

□ □

Type

DC Syringe filter
(2.8cm²)

□ □

Quantity/Package

N9 Only for syringe filter,
9pcs / pk
(3 batches, 3 pcs per batch)
N1 1pcs / pk

P

Market

P Biopharmaceutical



Viruclear™ VF Plus Removal Plastic-housing Cassette

V F

Application

VF Virus filtration

P

P Plus

E

Material of
Membrane

E mPES,
hydrophilic modified
PES membrane

S

Configuration
Code

S Standard

□ □

Type

LA Lab
(0.017/0.07m²)
FL Flow
(0.22/0.50/1.50m²)

□ □ □

Type

002 0.017m²
008 0.07m²
026 0.22m²
050 0.50m²
150 1.50m²

□ □

Quantity/
Package

N1 1pcs / pk

P

Market

P Biopharmaceutical



Ordering Information

Viruclear™ VF Plus X Virus Removal Syringe Filters with PVDF Shell

V F

Application

VF Virus filtration

X

X Plus X

E S

Material of Membrane

ES mPES, hydrophilic modified PES membrane

□ □

Type

DC Syringe filter (2.8cm²)

□ □

Quantity/Package

N9 Only for syringe filter, 9pcs / pk (3 batches, 3 pcs per batch)

N1 1pcs / pk

P

Market

P Biopharmaceutical



Viruclear™ VF Plus X Removal Plastic-housing Cassette

V F

Application

VF Virus filtration

X

X Plus X

E

Material of Membrane

E mPES, hydrophilic modified PES membrane

S

Configuration Code

S Standard

□ □

Type

LA Lab (0.017/0.07m²)

FL Flow (0.22/0.50/1.50m²)

□ □ □

Type

002 0.017m²

008 0.07m²

026 0.22m²

050 0.50m²

150 1.50m²

□ □

Quantity/Package


N1 1pcs / pk

P

Market

P Biopharmaceutical





**Filtration
Separation
Purification**

cobetter[®]
— filtration —

Please contact us for more information

Hangzhou Cobetter Filtration Equipment Co.,Ltd.

Sales Add Floor 6, Building 1, Cobetter Sci-Tech Innovation Center, No. 699 Mingxing Road,
Beigan Street, Xiaoshan District, Hangzhou City, Zhejiang Province, China

Factory Cobetter Park, Heshang New Material Industry Park,
Xiaoshan District, Hangzhou 311265, China

Tel +86 400-070-4266

Fax +86-571-87704256

www.cobetter.com

vipur

Purification Solutions

Your local distributor:

vipur GmbH

Mariahilfer Straße 113/18

1060 Vienna

AUSTRIA

+43 1 99 74 281

office@vipur.at

www.vipur.at